

EARTH ORTHO PVT LTD

Keeping You in Action

www.earthortho.in



info@earthortho.com




+91 - 72111 33303 Global Business
+91 - 72111 33305 Domestic Business



Corporate Office:


24-25, 3rd Floor, Shri Siddheshwar Plaza, Near Sakar Complex, New Vip Road,
Vadodara - 390022, Gujarat - India.



EARTH ORTHO PVT LTD

Keeping You in Action

LOCKING SYSTEM



2.7 MM / 2.7 MM VA - TWO COLUMN VOLAR DISTAL RADIUS PLATE



HEAD HOLE : 6, 7
SHAFT HOLES : 2, 3, 4, 5, 6, 7, 8
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS PLATE STRAIGHT - DORSAL



HEAD HOLE : -
SHAFT HOLES : 5, 6
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP "T" DISTAL RADIUS PLATE RIGHT ANGLE - DORSAL



HEAD HOLE : 3
SHAFT HOLES : 3, 4
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM CORTICAL

2.4 MM LCP "L" DISTAL RADIUS PLATE - DORSAL



HEAD HOLE : 2 -
SHAFT HOLES : 3, 4
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS PLATE - DORSAL



HEAD HOLE : 3
SHAFT HOLES : 3, 4
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP "L" DISTAL RADIUS PLATE OBLIQUE ANGLE - DORSAL



HEAD HOLE : 3
SHAFT HOLES : 3, 4
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS VOLAR PLATE - VOLAR



HEAD HOLE : 5
SHAFT HOLES : 3, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS EXTRA ARTICULAR PLATE - VOLAR



HEAD HOLE : 5
SHAFT HOLES : 3, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS EXTRA ARTICULAR PLATE - VOLAR



HEAD HOLE : 4
SHAFT HOLES : 3, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS BUTTRESS PLATE - VOLAR



HEAD HOLE : 5
SHAFT HOLES : 3
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS COLUMN VOLAR PLATE - VOLAR



HEAD HOLE : 8
SHAFT HOLES : 3, 4, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP DISTAL RADIUS COLUMN VOLAR PLATE - VOLAR



HEAD HOLE : 9
SHAFT HOLES : 3, 4, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP RADIAL HEAD NECK PLATE



HEAD HOLE : -
SHAFT HOLES : 2, 3, 4, 5
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL

2.4 MM LCP RADIAL HEAD RIM PLATE



HEAD HOLE : -
SHAFT HOLES : 2, 3, 4, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 MM LCP,
2.7 MM LCP, 2.7 MM CORTICAL



2.7 / 3.5 MM PERI - LOCK VOLAR DISTAL RADIUS LOCKING PLATE



HOLE : 3, 4, 5, 6, 7, 8
SHAFT HOLES : -
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
2.7 MM LCP, 3.5 MM CORTICAL

2.4 / 2.7 MM VA - LCP RIM DISTAL RADIUS PLATE



HEAD HOLE : 6, 7
SHAFT HOLES : 2, 3, 4, 5, 6, 7, 8
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : -

2.4 / 2.7 MM VA - DISTAL RADIUS META RIM PLATE



HOLE : 7
SHAFT HOLES : 3, 4, 5
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.4 / 2.7 MM LCP,
3.5 MM CORTICAL

5.0 MM TROCHANTRIC CABLE LCP PLATE



HOLE : -
SHAFT HOLES : 2, 4, 6, 8, 10, 12
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 5.0 MM LCP,
4.5 MM CORTICAL

2.4 MM LCP SCREW SELF TAPPING




DRILL BIT



LENGTH : 8 TO 40 MM (DIFF. 2 MM)
DRILL BIT : 1.8 MM
TYPE : SS / TT
DRILL SLEEVE : 1.8 / 2.4 MM

2.7 MM LCP SCREW SELF TAPPING




DRILL BIT



LENGTH : 8 TO 40 MM (DIFF. 2 MM)
DRILL BIT : 2.0 MM
TYPE : SS / TT
DRILL SLEEVE : 2.0 / 2.7 MM

2.7 MM CORTICAL SCREW SELF TAPPING



DRILL BIT



LENGTH : 8 TO 40 MM (DIFF. 2 MM)
DRILL BIT : 1.8 MM
TYPE : SS / TT
DRILL SLEEVE : 2.0 / 2.7 MM
BONE TAP : 2.7 MM

- 3.5 MM LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITHOUT LATERAL EXTENSION



HOLES : 6, 7, 8
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

- 3.5 MM LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION



HOLES : 3, 4, 5, 6, 7, 8
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 LHS, 3.5 LCP
 3.5 MM CORTICAL

- 3.5 MM LCP CLAVICLE HOOK PLATE



HOLES : 4, 5, 6, 7
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 DEPTH : 12, 15, 18
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

- 3.5 MM LCP PROXIMAL HUMERUS PLATE (SHORT PHILOS)




HOLES : 3, 4, 5
 ANGLED : -
 TYPE : SS / TT
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

- 3.5 MM LCP PROXIMAL HUMERUS PLATE (LONG PHILOS)




HOLES : 3, 4, 5, 6, 8, 10, 12
 ANGLED :-
 TYPE : SS / TT
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

- 3.5 MM LCP PERIARTICULAR PROXIMAL HUMERUS PLATE



HOLES : 2, 4, 6, 8, 10, 12, 14
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

- 3.5 MM LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE




HOLES : 4, 6, 8, 10, 12, 14
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

- 2.7 / 3.5 MM LCP POSTEROLATERAL DISTAL HUMERUS PLATE WITH SUPPORT



HOLES : 3, 5, 7, 9, 11, 13
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 MM LHS,
 3.5 MM LCP, 3.5 MM CORTICAL

- 2.7 / 3.5 MM LCP POSTEROLATERAL DISTAL HUMERUS PLATE WITHOUT SUPPORT



HOLES : 3, 5, 7, 9, 11, 13
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 LHS, 3.5 LCP
 3.5 MM CORTICAL

- 3.5 MM LCP MEDIAL DISTAL HUMERUS PLATE



HOLES : 3, 5, 7, 9, 11, 13
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 MM LHS,
 3.5 MM LCP, 3.5 MM CORTICAL

- 3.5 MM LCP METAPHYSEAL PLATE FOR DISTAL MEDIAL HUMERUS




HOLES : 7, 9, 11, 13
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 LHS, 3.5 LCP
 3.5 MM CORTICAL

- 3.5 MM LCP OLECRANON PLATE



HOLES : 4, 6, 8, 10, 12
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 LHS, 3.5 LCP
 3.5 MM CORTICAL

- 3.5 MM LCP PROXIMAL HUMERUS PLATE (LPHP)



HOLES : 5, 6, 7, 8
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 MM LCP,
 3.5 MM LCP, 3.5 MM CORTICAL

- 3.5 MM LCP 1/3 TUBULAR PLATE/
3.5 MM COMBI HOLE RECON PLATE




HOLES : 3 TO 12
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 MM LCP,
 3.5 MM LCP, 3.5 MM CORTICAL

- 3.5 MM LCP RECONSTRUCTION PLATE (STRAIGHT)



HOLES : 4 TO 16
 ANGLED : -
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 MM LCP,
 3.5 MM LCP, 3.5 MM CORTICAL

- 3.5 MM LCP METAPHYSEAL PLATE




HOLES : 3, 4, 5, 6, 7, 8, 9, 11, 13
 ANGLED : -
 TYPE : SS / TT
 SUITABLE SCREW : 2.7 MM LCP,
 3.5 MM LCP, 3.5 MM CORTICAL

● **3.5 MM LCP POSTERIOR MEDIAL PROXIMAL TIBIA PLATE**



HOLES : 1, 2, 4, 6, 8, 10
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP MEDIAL PROXIMAL TIBIA PLATE**




HOLES : 4, 6, 8, 10, 12, 14, 16
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP PROXIMAL TIBIA PLATE (RAFT)**




HOLES : 4, 6, 8, 10, 12, 14, 16
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP LOW BEND MEDIAL DISTAL TIBIA PLATE**



HOLES : 4, 6, 8, 10, 12, 14
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP MEDIAL DISTAL TIBIA PLATE WITH TAB**



HOLES : 4, 6, 8, 10, 12, 14
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP ANTEROLATERAL DISTAL PLATE**




HOLES : 5, 7, 9, 11, 13, 15, 17
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP T - PLATE RIGHT ANGLE**



HOLES : 3, 4, 5, 6, 7, 8
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP T - PLATE OBLIQUE ANGLE**




HOLES : 3, 4, 5, 6, 7, 8
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP CLOVERLEAF PLATE**




HOLES : 4, 5, 6, 7, 8, 9, 10, 11, 12
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP CALCANEAL PLATE**




HOLES : 64, 69, 76, 81
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **2.7 / 3.5 MM LCP LATERAL DISTAL FIBULA PLATE**




HOLES : 3, 4, 5, 6, 7, 8, 9, 10, 11
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.7 MM LCP,
 3.5 MM LCP, 3.5 MM CORTICAL

● **2.7 / 3.5 MM LCP POSTEROLATERAL DISTAL FIBULA PLATE**




HOLES : 3, 4, 5, 6, 7, 8, 9, 10, 11
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 2.7 MM LHS,
 3.5 MM LCP, 3.5 MM CORTICAL

● **3.5 MM LCP MEDIAL ANTERIOR CLAVICLE PLATE**




HOLES : 6, 7, 8
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP PILON PLATE**



HOLES : 7, 9
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP SMALL PLATE**



HOLES : 4, 5, 6, 7, 8, 9, 10, 11,
 12, 13, 14
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

● **3.5 MM LCP SMALL HOOK PLATE**



HOLES : 3, 4, 5, 6, 7, 8
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL

3.5 MM LOCKING SYSTEM
SUPPORTED SCREWS & DRILL BITS



3.5 MM LCP SCREW SELF TAPPING



LENGTH : 10 TO 60 MM (DIFF. 2 MM)
LENGTH : 65 TO 70 MM (DIFF. 5 MM)
TYPE : SS / TT
DRILL BIT : 2.8 MM

2.7 MM LCP SCREW SELF TAPPING



LENGTH : 8 TO 40 MM (DIFF. 2 MM)
DRILL BIT : 2.0 MM
TYPE : SS / TT
DRILL SLEEVE : 2.0 / 2.7 MM

3.5 MM CORTICAL SCREW SELF TAPPING



LENGTH : 10 TO 70 MM (DIFF. 2 MM)
DRILL BIT : 2.0 MM
TYPE : SS / TT
DRILL SLEEVE : 2.5 / 3.5 MM
BONE TAP : 3.5 MM

4.0 MM CANCELLOUS SCREW FULLY
THREADED




LENGTH : 10 TO 30 MM (DIFF. 2 MM)
LENGTH : 35 TO 75 MM (DIFF. 5 MM)
TYPE : SS / TT
DRILL & TAP SLEEVE : 2.0 / 4.0 MM
BONE TAP : 4.0 MM

4.0 MM CANCELLOUS SCREW PARTIALLY
THREADED



LENGTH : 10 TO 30 MM (DIFF. 2 MM)
LENGTH : 35 TO 75 MM (DIFF. 5 MM)
TYPE : SS / TT
DRILL & TAP SLEEVE : 2.0 / 4.0 MM
BONE TAP : 4.0 MM

3.5 MM LCP CANCELLOUS SCREW FULLY
THREADED



LENGTH : 10 TO 30 MM (DIFF. 2 MM)
LENGTH : 35 TO 70 MM (DIFF. 5 MM)
TYPE : SS / TT
STEP DRILL BIT : 2.5 / 2.8 MM
DRILL SLEEVE : 2.5 / 2.8 / 3.5 MM

3.5 MM MALLEOLAR SCREW SELF TAPPING



LENGTH : 20 TO 70 MM (DIFF. 5 MM)
TYPE : SS / TT



2.7 / 3.5 MM LCP INSTRUMENT
SYSTEM



SMALL FRAGMENT
INSTRUMENT SET
(BASIC + LCP)




For more information
+91 72111 33303 (Export)
+91 72111 33305 (Domestic)

ORDER NOW

www.earthortho.in

● 4.5 / 5.0 MM LCP DISTAL FEMUR PLATE



HOLES : 5, 7, 9, 11, 13, 15
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 / 6.5 MM LCP PROXIMAL FEMUR PLATE



HOLES : 4, 6, 8, 10, 12, 14
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 / 6.5 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 / 6.5 MM LCP PROXIMAL FEMUR HOOK PLATE



HOLES : 4, 6, 8, 10, 12, 14, 16
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 / 6.5 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP LATERAL DISTAL FEMUR PLATE



HOLES : 4
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP LATERAL HIGH TIBIA PLATE




HOLES : 3
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP PROXIMAL LATERAL TIBIA PLATE



HOLES : 5, 7, 9, 11, 13, 15
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP PROXIMAL TIBIA PLATE



HOLES : 4, 6, 8, 10, 12, 14, 16
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP NARROW PLATE




HOLES : 4, 6, 8, 10, 12, 14, 16
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP BROAD PLATE




HOLES : 5, 6, 7, 8, 9, 10, 11, 12,
 13, 14, 15, 16, 18
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP RECONSTRUCTION PLATE (STRAIGHT)



HOLES : 4, 5, 6, 7, 8, 9, 10, 11, 12,
 13, 14, 15, 16
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP METAPHYSEAL PLATE




HOLES : 3, 4, 5, 6, 7, 8, 9, 10, 11, 13
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP T - BUTTRESS PLATE



HOLES : 4, 5, 6, 7, 8, 9, 10, 11, 12
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP L - BUTTRESS PLATE




HOLES : 4, 5, 6, 7, 8, 9, 10, 11, 12
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP MEDIAL HIGH TIBIA PLATE



HOLES : 4
ANGLED : -
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 / 6.5 MM LCP CONDYLAR BUTTRESS PLATE



HOLES : 6, 8, 10, 12, 14
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 / 6.5 MM CANCELLOUS

● 4.5 / 5.0 LCP T - PLATE



HOLES : 4, 5, 6, 7, 8, 9, 11, 12
ANGLED : LEFT / RIGHT
TYPE : SS / TT
SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM LCP LATERAL DISTAL FEMUR PLATE



HOLES : 5, 7, 9, 11, 13, 15
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 4.5 MM CORTICAL,
 5.0 MM LCP, 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM 130° & 135° LCP DHS PLATE LONG BARREL




HOLES : 3, 4, 5, 6, 7, 8, 9, 10, 11, 12
 ANGLED : -
 TYPE : SS / TT
 SUITABLE SCREW : 4.5 / 5.0 MM,
 12.5 MM DHS CANCELLOUS

● 4.5 / 5.0 MM PROXIMAL FEMUR LCP PLATE



HOLES : 4, 6, 8, 10, 12, 14
 ANGLED : LEFT / RIGHT
 TYPE : SS / TT
 SUITABLE SCREW : 5.0 MM LCP,
 5.0 MM CANCELLOUS

● 4.5 / 5.0 MM 130° & 135° LCP DHS PLATE SHORT BARREL



HOLES : 3, 4, 5, 6, 7, 8, 9, 10, 11, 12
 ANGLED : -
 TYPE : SS / TT
 SUITABLE SCREW : 4.5 / 5.0 MM LCP,
 12.5 MM DHS SCREW

● 4.5 / 5.0 / 6.5 MM 95° LCP LCP DCS PLATE



HOLES : 4, 5, 6, 7, 8, 9, 10, 11, 12, 14
 ANGLED : -
 TYPE : SS / TT
 SUITABLE SCREW : 4.5 / 5.0 MM LCP,
 6.5 MM CANCELLOUS

● 12.5 MM DHS / DCS COMPRESSION LAG SCREW




LENGTH : 40 TO 120 MM (DIFF. 5 MM)
 TYPE : SS / TT

● 4.5 / 5.0 MM TROCHANTER STABILIZING PLATE FOR DHS - ADJUSTABLE




SUITABLE SCREW : 4.5 / 5.0 MM LCP,
 12.5 MM DHS SCREW

● 3.5 MM LCP WRIST FUSION PLATE



STANDARD BEND : 8 HOLES
 SHORT BEND : 8 HOLES
 STRAIGHT BEND : 9 HOLES
 SUITABLE SCREW : 3.5 MM LCP,
 3.5 MM CORTICAL


● 5.0 MM LCP SCREW SELF TAPPING



DRILL BIT 

LENGTH : 14 TO 60 MM (DIFF. 2 MM)
 LENGTH : 65 TO 100 MM (DIFF. 5 MM)
 TYPE : SS / TT
 DRILL BIT : 4.3 MM
 DRILL SLEEVE : 4.3 / 5.0 MM

● 4.5 MM CORTICAL SCREW SELF TAPPING



DRILL BIT 

LENGTH : 14 TO 70 MM (DIFF. 2 MM)
 TYPE : SS / TT
 DRILL BIT : 3.2 MM
 DRILL & TAP SLEEVE : 3.2 / 4.5 MM
 BONE TAP : 4.5 MM

● 6.5 MM CANCELLOUS SCREW 16 MM THREADED



DRILL BIT 

LENGTH : 25 TO 110 MM (DIFF. 5 MM)
 TYPE : SS / TT
 DRILL BIT : 3.2 MM, 4.5 MM
 DRILL & TAP SLEEVE : 3.2 / 6.5 MM
 BONE TAP : 6.5 MM

● 6.5 MM CANCELLOUS SCREW 32 MM THREADED



DRILL BIT 

LENGTH : 40 TO 110 MM (DIFF. 5 MM)
 TYPE : SS / TT
 DRILL BIT : 3.2 MM, 4.5 MM
 DRILL & TAP SLEEVE : 3.2 / 6.5 MM
 BONE TAP : 6.5 MM

● 6.5 MM CANCELLOUS SCREW FULLY THREADED



DRILL BIT 

LENGTH : 25 TO 110 MM (DIFF. 5 MM)
 TYPE : SS / TT
 DRILL BIT : 3.2 MM, 4.5 MM
 DRILL & TAP SLEEVE : 3.2 / 6.5 MM
 BONE TAP : 6.5 MM


● 5.0 MM LCP CANCELLOUS SCREW FULLY THREADED



DRILL BIT 

LENGTH : 25 TO 110 (DIFF. 5 MM)
 TYPE : SS / TT
 STEP DRILL BIT : 3.5 MM, 4.3 MM
 DRILL SLEEVE : 3.5 / 6.5 MM

● 6.5 MM CANCELLOUS SCREW FULLY THREADED



LENGTH : 25 TO 110 (DIFF. 5 MM)
 TYPE : SS / TT
 STEP DRILL BIT : 3.5 MM, 4.3 MM
 DRILL SLEEVE : 3.5 / 6.5 MM

● 4.5 MM MALLEOLAR SELF TAPPING



LENGTH : 20 TO 100 MM (DIFF. 5 MM)
 TYPE : SS / TT

EARTH ORTHO PVT LTD
Keeping You in Action

**LARGE FRAGMENT
INSTRUMENT SET
(BASIC + LCP)**

Tray 1

Tray 2

Tray 3

Screw Box


For more information
+91 72111 33303 (Export)
+91 72111 33305 (Domestic)

ORDER NOW

www.earthortho.in



IMPLANT TRAY



SCREW TRAY

