




EARTH ORTHO PVT LTD

Keeping You in Action



INTRAMEDULLARY NAILS

● HUMERUS NAIL CANNULATED

Ø6.5 MM - 3.4 MM Interlocking Bolt

Ø7.5 MM - 3.4 MM Interlocking Bolt

Distal Holes

3 Distal Holes



3 Proximal Holes

Proximal Holes

Ø6.5 MM - 3.4 MM Interlocking Bolt

Ø7.5 MM - 3.4 MM Interlocking Bolt

LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS / TT
DIA Ø	6.5, 7.3 MM

● TIBIA NAIL CANNULATED

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 4.9 MM Interlocking Bolt

Distal Holes

3 Distal Holes



3 Proximal Holes

Proximal Holes

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 4.9 MM Interlocking Bolt

LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0 MM

● EXPERT TIBIA NAIL CANNULATED

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 MM to 11.0 MM - 3.9 & 4.9 MM Interlocking Bolt

Distal Holes

4 Distal Holes



5 Proximal Holes

Proximal Holes

Ø8.0 MM - 3.9 MM Interlocking Bolt

Ø9.0 to 11.0 MM - 4.9 MM Interlocking Bolt

LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0 MM

Intramedullary Nails

● FEMUR NAIL CANNULATED


Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt



LENGTH	320, 340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM

● 130° & 135° SHORT PROXIMAL FEMUR NAIL CANNULATED

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt



LENGTH	250 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM


● 130° & 135° TROCHANTERIC PROXIMAL FEMUR NAIL CANNULATED

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt




LENGTH	180 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM

● 130° SHORT PROXIMAL FEMUR NAIL ANTIROTATION - II CANNULATED



● 130° TROCHANTERIC PROXIMAL FEMUR NAIL ANTIROTATION - II CANNULATED




● 130° & 135° LONG PROXIMAL FEMUR NAIL CANNULATED

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt

Distal Holes



2 Distal Holes



2 Proximal Holes



Proximal Holes

Ø9.0 MM TO 12.0 MM - 6.4 MM CANNULATED PFN HIP SELF TAPPING
Ø9.0 MM TO 12.0 MM - 8.0 MM CANNULATED PFN FEMORAL NECK
SCREW SELF TAPPING

LENGTH	340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	LEFT / RIGHT


● 130° LONG PROXIMAL FEMUR NAIL ANTIROTATION - II CANNULATED

Ø9.0 MM to 12.0 MM - 4.9 MM Interlocking Bolt

Distal Holes



2 Distal Holes



1 Proximal Holes



Proximal Holes

Ø9.0 MM to 12.0 MM - 12.0 MM CANNULATED PFNA
- II BLADE


LENGTH	340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	LEFT / RIGHT

● SUPRA CONDYLAR DISTAL FEMUR NAIL CANNULATED



LENGTH	150, 200, 250, 300, 350, 400 MM
TYPE	SS / TT
DIA Ø	9.0, 10.0, 11.0, 12.0 MM
ANGLED	UNIVERSAL

● ELASTIC NAIL (TENS)



1.5 MM
2.0 MM
2.5 MM
3.0 MM
3.5 MM
4.0 MM
4.5 MM

LENGTH	440 MM
TYPE	SS / TT
DIA Ø	1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5 MM
ANGLED	UNIVERSAL

● ENDER'S NAIL



LENGTH	200, 220, 240, 260, 280, 300, 320, 340, 360, 380, 400, 420, 440 MM
TYPE	SS / TT
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0, 4.5 MM

● 12.0 MM CANNULATED PFNA - II BLADE



LENGTH	50 MM TO 120 MM (DIFF. 5 MM)
TYPE	SS / TT
DIA Ø	12.0 MM

● SQUARE NAIL FOR ULNA




LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS / TT
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0 MM

● SQUARE NAIL FOR RADIUS



LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS / TT
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0 MM

● RUSH NAIL



LENGTH	200, 220, 240, 260, 280, 300 MM
TYPE	SS / TT
DIA Ø	2.0, 2.5, 3.0, 3.5, 4.0 MM

● K - NAIL



LENGTH	280, 300, 320, 340, 360, 380 MM
TYPE	SS / TT
DIA Ø	8.0, 9.0, 10.0, 11.0, 12.0, 13.0 MM

● END CAP - TIBIA / FEMUR / PEN / SUPRA



LENGTH	0, 5, 10, 15 MM
TYPE	SS / TT

● 3.4 MM INTERLOCKING BOLT SELF TAPPING



LENGTH	14 TO 50 MM (DIFF. 2 MM)
TYPE	SS / TT
DRILL BIT	2.7 MM

● 3.9 MM INTERLOCKING BOLT SELF TAPPING



LENGTH	14 TO 80 MM (DIFF. 2 MM)
TYPE	SS / TT
DRILL BIT	3.2 MM

● 4.9 MM INTERLOCKING BOLT SELF TAPPING



LENGTH	18 TO 80 MM (DIFF. 5 MM)
TYPE	SS / TT
DRILL BIT	4.0 MM

● 6.4 MM CANNULATED PFN HIP PIN SELF TAPPING



LENGTH	40 TO 110 MM (DIFF. 5 MM)
TYPE	SS / TT
REAMER	6.4 MM
GUIDE WIRE	2.5 MM

● 8.0 MM CANNULATED PFN FEMORAL NECK SCREW SELF TAPPING



LENGTH	40 TO 110 MM (DIFF. 5 MM)
TYPE	SS / TT
REAMER	8.0 MM
GUIDE WIRE	2.5 MM

PFN II INSTRUMENT SET



   FOLLOW US



PFN II INSTRUMENT SET

ORDER NOW



+91 72111 33303

+91 72111 33305

 www.earthortho.in

 info@earthortho.com





TENS NAIL INSTRUMENT



EARTH ORTHO PVT LTD
Keeping You in Action



TENS NAIL **INSTRUMENT SYSTEM**



More Info Call Us:



+91 72111 33303

+91 72111 33305

www.earthortho.in



TENS NAIL INSTRUMENT



EARTH ORTHO PVT LTD
Keeping You in Action

   FOLLOW US

 **CONTACT US**

+91 72111 33303

(Export)

+91 72111 33305

(Domestic)

**ORDER
NOW**

"WE ARE COMMITTED
TO **CUSTOMER SATISFACTION**"



info@earthortho.com



www.earthortho.com



**Tibia Femur
Instrument set**





EARTH ORTHO PVT LTD

 **Keeping You in Action**



info@earthortho.com



www.earthortho.in



+91 - 72111 33303 Global Business
+91 - 72111 33305 Domestic Business



Corporate Office:

**24-25, 3rd Floor, Shri Siddheshwar Plaza, Near Sakar Complex, New Vip Road,
Vadodara - 390022, Gujarat - India.**