




JOINT REPLACEMENT

- **BIPOLAR PROSTHESIS
(FENESTRATED)**

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **BIPOLAR PROSTHESIS
(NON-FENESTRATED)**

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- AUSTIN MOORE HIP PROSTHESIS**

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- AUSTIN MOORE HIP PROSTHESIS (NARROW STEM)**

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **AUSTIN MOORE HIP PROSTHESIS**
(LONG STEM)

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **THOMPSON HIP PROSTHESIS**

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **THOMPSON HIP PROSTHESIS
STANDARD**
(LONG STEM)

SIZE	37 MM TO 55 MM (DIFF. 1 MM)
MATERIAL	SS



- **SHOULDER PROSTHESIS**

SIZE	31 MM TO 41 MM (DIFF. 2 MM)
MATERIAL	SS
DIA Ø	6.0, 0.7, 8.0 MM



● CEMENTED STEM (POLISHED)

MATERIAL	SIZE	LENGTH
SS 316L	NO - 8	115 MM
	NO - 9	130 MM
	NO - 10	140 MM
	NO - 11	145 MM
	NO - 12	150 MM
	NO - 13	155 MM
	NO - 14	160 MM
	NO - 15	165 MM

(12 / 14 TAPPER & NECK ANGLE 135°)
(SELF CENTERING)



● UNCEMENTED STEM (H. A. COATED)



MATERIAL	SIZE	LENGTH
TATINUM	NO - 8	115 MM
	NO - 9	130 MM
	NO - 10	140 MM
	NO - 11	145 MM
	NO - 12	150 MM
	NO - 13	155 MM
	NO - 14	160 MM
	NO - 15	165 MM

(12 / 14 TAPPER & NECK ANGLE 135°)
(SELF CENTERING)


● UNCEMENTED STEM (POROUS COATED)

MATERIAL	SIZE	LENGTH
SS 316L	NO - 8	115 MM
	NO - 9	130 MM
	NO - 10	140 MM
	NO - 11	145 MM
	NO - 12	150 MM
	NO - 13	155 MM
	NO - 14	160 MM
	NO - 15	165 MM

(12 / 14 TAPPER & NECK ANGLE 135°)
(SELF CENTERING)



● MODULAR BIPOLAR CUP



SIZE	38 MM TO 56 MM
MATERIAL	SS 316L

(SELF CENTERING)

● INNER HEAD (POLISHED)

MATERIAL	DIA - Ø 28 MM (12 / 14 TAPPER)			
SS 316L	-3.5	0.0	+3.5	+7.0



● CENTRELIZER


● CEMENT RETRACTOR



- **UNCEMENTED ACETABULAR CUP
WITH LINER UHMWPW 10° & 180°
(H. A. COATED)**

SIZE	44 MM, 46 MM, 48 MM, 50 MM, 52 MM, 54 MM, 56 MM, 58 MM, 60 MM, 62 MM, 64 MM
MATERIAL	TITANIUM (H.A. COATED)

(SELF CENTERING)



- **UNCEMENTED ACETABULAR CUP
WITH LINER UHMWPW 10° & 180°
(POROUS COATED)**



SIZE	44 MM, 46 MM, 48 MM, 50 MM, 52 MM, 54 MM, 56 MM, 58 MM, 60 MM, 62 MM, 64 MM
MATERIAL	TITANIUM (POROUS COATED)

(SELF CENTERING)

● CEMENTED ACETABULAR CUP 10° & 180°

SIZE	44 MM, 46 MM, 48 MM, 50 MM, 52 MM, 54 MM, 56 MM, 58 MM, 60 MM, 62 MM, 64 MM
MATERIAL	UHMWPW

(SELF CENTERING)







● ACETABULAR CUP SCREW

SIZE Ø6.5 MM	20 MM TO 50 MM (DIFF. 5 MM)
MATERIAL	TITANIUM



INSTRUMENT SET



INSTRUMENT SET







HIP REPLACEMENT INSTRUMENT SET

Contact Us

+91 72111 33303

+91 72111 33305



www.earthortho.in



info@earthortho.com



INSTRUMENT SET

   FOLLOW US



BIOPOLARE INSTRUMENT SET



 ORDER NOW

 +91 72111 33303
+91 72111 33305

 info@earthortho.com  www.earthortho.in



INSTRUMENT SET



EARTH ORTHO PVT LTD
Keeping You in Action

   FOLLOW US



THOMPSON INSTRUMENT SET

 +91 72111 33303, +91 72111 33305


www.earthortho.in



   FOLLOW US



SHOULDER PROSTHESIS INSTRUMENT SET



JOIN US

 info@earthortho.com

 www.earthortho.in




+91 72111 33303

+91 72111 33305



TOTAL HIP REPLACEMENT





EARTH ORTHO PVT LTD

 **Keeping You in Action**



info@earthortho.com



www.earthortho.in



+91 - 72111 33303 Global Business
+91 - 72111 33305 Domestic Business



Corporate Office:

**24-25, 3rd Floor, Shri Siddheshwar Plaza, Near Sakar Complex, New Vip Road,
Vadodara - 390022, Gujarat - India.**